

# Configura® Lite 2

## The entry-level chair with tilt-in-space

The Lite 2 is the most cost-effective chair in the **Configura** range, featuring a single motor tilt-in-space mechanism, height and depth adjustment and a high-risk pressure cushion as standard.



|  | Small           | Medium          | Large |
|--|-----------------|-----------------|-------|
| Overall chair width                        | 750 mm (29.5")  | 800 mm (31.5")  | n/a   |
| Overall height                             | 1080 mm (42.5") | 1115 mm (43.9") | n/a   |
| Overall length                             | 885 mm (34.8")  | 940 mm (37")    | n/a   |
| Overall chair depth with legs elevated     | 1400 mm (55")   | 1400 mm (55")   | n/a   |
| Internal seat width                        | 457 mm (18")    | 508 mm (20")    | n/a   |
| Seat height (from footplate if applicable) | n/a             | n/a             | n/a   |
| Floor to seat height                       | 406 mm (16")    | 457 mm (18")    | n/a   |
| Internal seat depth Position 1             | 457 mm (18")    | 457 mm (18")    | n/a   |
| Internal seat depth Position 2             | 508 mm (20")    | 508 mm (20")    | n/a   |
| Backrest height from seat                  | 608 mm (27")    | 608 mm (27")    | n/a   |
| Chair to wall/furniture                    | 350 mm (13.8")  | 350 mm (13.8")  | n/a   |
| Backrest angle – std                       | 66°             | 66°             | n/a   |
| Backrest angle – recline                   | 48°             | 48°             | n/a   |
| Rise angle – std                           | 23°             | 23°             | n/a   |
| Rise angle – med rise                      | n/a             | n/a             | n/a   |
| Rise angle – flat                          | n/a             | n/a             | n/a   |
| Rise height at front of seat – std         | 510 mm (20.1")  | 545 mm (21.5")  | n/a   |
| Rise height at front of seat – med         | 570 mm (22.4")  | 620 mm (24.4")  | n/a   |
| Rise height at front of seat – flat        | 660 mm (26")    | 710 mm (28")    | n/a   |
| TIS angle                                  | 26°             | 26°             | n/a   |
| Standard seat angle in neutral             | 3°              | 3°              | n/a   |
| Forward tilt angle                         | n/a             | n/a             | n/a   |
| Armrest height                             | 180 mm (7.1")   | 180 mm (7.1")   | n/a   |
| Weight capacity                            | 125 kg (20 st)  | 160 kg (25 st)  | n/a   |
| Leg rest weight capacity                   |                 |                 | n/a   |
| Overall weight                             | 65 kg (143 lbs) | 68 kg (150 lbs) | n/a   |