

# ConfiguraSolus®

## The portable riser recliner chair compatible with all transfer methods

The Solus can be used as a riser recliner and a porter chair. It can be upgraded as a user's needs change to ensure the chair is always suitable for the user. The adjustability also means it can be easily reused by another user.



	Min	Max
Seat height	400 mm (16")	550 mm (22")
Seat width	400 mm (16")	550 mm (22")
Seat depth	400 mm (16")	550 mm (22")
Armrest height	95 mm (3.7")	95 mm (3.7")
Overall width	700 mm (27.5")	770 mm (30.3")
Overall height - upright	1100 mm (43.3")	1240 mm (48.8")
Overall length - upright	880 mm (34.6")	930 mm (36.6")
Rise height to front of chair (standard rise)	540 mm (21.2")	680 mm (26.7")
Rise height to front of chair (medium rise)	600 mm (23.6")	740 mm (29.1")
Total length fully reclined	1550 mm (61")	1720 mm (67.7")
Chair to wall/furniture	600 mm (23.6")	770 mm (30.3")
Standard seat angle	6°	6°
Tilt-in-space seat angle	28°	28°
Rise angle - standard position	18°	18°
Rise angle - medium position	12°	12°
Standard backrest angle	95°	95°
Recline backrest angle	125°	125°