

# Accora

At Accora, we create and develop assistive products to help people lead as full and safe a life as possible. We aim to safeguard and enrich people's lives, bringing dignity and creating possibility.

## FloorBeds



### EMPRESA®

- Lifting range: 3.9 to 31.5 inches
- BodyMove technology supporting pressure ulcer prevention
- Option to extend the width and length by 4 inches
- Available in 100+ design options
- Weight capacity: 600 lbs

NSB-2-FL6-000US



### CONTESA®

- Lifting range: 3.9 to 29 inches
- Full range of profiling options for care, comfort, and positioning
- 36- or 42-inch width with an optional 4-inch length extension
- Available in 100+ design options
- Weight capacity: 500 lbs

36" - NSB-0-FL8-000US

42" - NSB-W-0-FL8-000US



### FloorBed® 1

- Lifting range: 2.8 to 25.6 inches
- Full range of profiling functions to aid care and positioning
- Weight capacity: 408 lbs
- 42-inch width option with 575 lbs weight capacity available

36" - NSB-0-FL1-200US

42" - NSB-0-FL5-200US

## Seating




### ConfiguraAdvance®

With pressure-relieving functions and infinite adjustability, easily adjust seat depth and width, footplate height and angle and backrest cushions and angle, all without tools.

CHAIR-0-SC3-050

# Therapy surfaces

## Allevia® Comfort




### Foam

- Die-cut high-density foam promotes consistent pressure distribution
- Firm edges
- Sloped memory foam heel section
- Multi-directional stretch nylon cover
- Weight capacity: 400 lbs

36" x 80" – ACC400-0-FM1-100    36" x 84" – ACC400-0-FM1-200  
42" x 80" – ACC400-0-FM2-100    42" x 84" – ACC400-0-FM2-200

## Allevia® Prima




### Foam

- Gel-infused memory foam top layer
- Die cut medium-density foam middle layer with a high-density foam base layer
- Firm edges
- Sloped memory foam heel section
- Multi-directional stretch nylon cover
- Weight capacity: 500 lbs

36" x 80" – ACC600-0-FM1-100    36" x 84" – ACC600-0-FM1-200  
42" x 80" – ACC600-0-FM2-100    42" x 84" – ACC600-0-FM2-200

## Allevia® Duo



### Hybrid

- Ten connected cells redistribute pressure as the user moves
- Can be upgraded to alternating powered therapy with an optional pump
- Gel-infused memory foam top with a high-density foam middle layer and sloped heel section
- Four-way stretch nylon cover
- Weight capacity: 36" – 500 lbs    42" – 600 lbs

**Mattress – ACC800-0-AF1-100**

**Mattress with pump – ACC800P-0-AF1-100**

## Allevia® Sense



### Air

- 7-inch air cells that can be deflated or repositioned
- 4.5-inch heel slope cells that can be deflated or repositioned
- Auto-Sense™ technology automatically detects patient weight
- Functions: Alternating or constant low-pressure therapy, seat inflation, and AutoFirm mode for care tasks
- Cell-in-cell design maintains support during power loss or transport
- Ultra-rapid CPR valve enables quick emergency deflation
- Weight capacity: 550 lbs

36" – AIRMATP-0-AM4-000

42" – AIRMATP-0-AM4-100