

The most versatile riser recliner

Building on the Lite 2, the Comfort is a dual motor tilt-in-space chair with independent no-shear backrest recline, offering standard and flat rise options. The subframe is also designed to be compatible with a range of transfer aids.



Medium



Key features:



Flat rise (option)



Tilt-in-space with legrest elevation



Backrest recline



Adjustable seat depth



Braked castors



Integrated self-supply overlay

Aperture dimensions
17.5" x 17.5" x 3"



Pressure-reducing fabric

As standard on backrest and armrest.




Compatible with transfer aids

Cricket, ReTurn 7500 and Molift Raiser.



'Reduced shear' recline

Maximises pelvic stability.



Depth adjustment options using backrest brackets

Reduces the need for depth adjustment kits.



Drop-down armrest

Available as an optional extra.



Vinyl fabric

Available as an optional extra.



Cocoon backrest

Available as an optional extra.

Chair size	Part number	Seat depth	Seat height	Seat width	Weight capacity
Small	CHAIR-0-CMI-140	450 mm (18")	375 mm (15")*	450 mm (18")	25 stone
			400 mm (16")		
		500 mm (20")	425 mm (17")*		
Medium	CHAIR-0-CMI-240	450 mm (18")	425 mm (17")*	500 mm (20")	25 stone
			450 mm (18")**		
		500 mm (20")	475 mm (19")*		
			500 mm (20")**		
Large	CHAIR-0-CMI-340	450 mm (18")	425 mm (17")*	550 mm (22")	25 stone
			450 mm (18")**		
		500 mm (20")	475 mm (19")*		
			500 mm (20")**		

** Includes 2" Height Adjustment Kit * Includes 1" Height Adjustment Kit ° The castors are removed to make the chair lower
 ■ Standard size ■ Seat depth is adjustable by 50 mm (2") as standard without extra parts