

# Community Bed



# CommunityBed

CommunityBed allows equipment services to combine two beds into one with its lifting range of 220 mm – 800 mm. By combining the lifting height of a conventional profiling bed (800 mm) and the lowness of a conventional low bed (220 mm), CommunityBed allows equipment service to simplify their range and reduce stock holding.

CommunityBed is fully compliant with BSEN 60601-2-52 and compatible with a broad range of accessories including integrated, padded mesh side rails.



## KEY FEATURES

- Lifting range of 220 mm – 800 mm
- All parts under 25 kg for ease of delivery
- Auto-regression on backrest
- With 6" mattress, side rails are 290 mm above top of mattress
- With 8" mattress, side rails are 240 mm above top of mattress which is above the 220 mm measurement required by BSEN 60601-2-52 standards

## TECHNICAL SPECIFICATION

|                      |  |
|----------------------|--|
| Dimensions:          | 930 x 2250 mm (36" x 89")  |
| Platform size:       | 900 x 2000 mm (36" x 79")  |
| Safe working load:   | 215 kg (34 stone)  |
| Max patient weight:  | 180 kg (28 stone)  |
| Lifting range:       | 220 - 800 mm (8.6" - 31.5")<br><i>(excl. mattress)</i>   |
| Profiling functions: | Backrest with Auto-regression, Kneebreak. Anti-Trendelenburg and Trendelenburg available as an option. |
| Part number:         | NSB-0-CB1-000  |



Easy to use handset



Excellent height for carers



Fully profiling



Bed lever




Integrated side rails - raised




Integrated side rails - lowered

# Accora Beds Comparison

## COMMUNITYBED



## FLOORBED®



### Injury risks when falling from bed

