

# ConfiguraAdvance<sup>®</sup>


**Comprehensive posture management**  
designed for high need users

Accora




With integrated size adjustments in every dimension and backrest upgrade options, the **Advance** is our most configurable and recyclable chair yet, making it the optimal choice for clients with complex and postural requirements.





## Size

Change the seat depth, width and height as well as the footplate height and angle according to the individual needs of the service user, all without tools.



## Pressure

Integrated and adaptable pressure management for your peace of mind.



### Height adjustment

Adjust seat height from 330 mm (13") to 530 mm (21")

### Depth adjustment

Adjust seat depth from 400 mm (16") to 550 mm (22")



### Width adjustment

Adjust seat width from 400 mm (16") to 550 mm (22").



### Integrated self-supply overlay

Fit your own pressure cushion.

### Pressure-reducing fabric

Continuous overlay for thigh support.



### BodyMove, reduced shear backrest recline

Maximises pelvic stability and assists in preventing tissue damage. Full length compensation on leg rest giving superior leg rest elevation and support.




## Posture

Infinitely configurable support options suitable for even complex posture management needs.

## Optional extras to complement the ConfiguraAdvance®



**Waterfall** backrest



**Lateral** supports



**Cocoon** backrest



**Forward tilt**



**Electric operation**



**Dropdown armrests**



**Removable** footplate

Features angle and height adjustment options as standard.



**Postural backrest upgrade**

Available as an option to support more complex postural requirements.



**Four point pelvic positioning belt**



**Armrest height**




**Cocoon backrest**



**Mattress infill**



**Junior Kit**



**External lateral supports**

Technical specification

Part number	Adjustment	Range	Weight capacity
CHAIR-0-SCI-030	Depth	400 mm (16") – 550 mm (22")	25 stone
	Height	330 mm (13") – 530 mm (21")	
	Width	350 mm (14") – 550 mm (22")	

**UK:** Accora Ltd, Charter House, Barrington Road, Orwell, Cambridge SG8 5QP  
**T:** +44 (0)1223 206 100 **E:** info@accora.care **W:** www.accora.care

**IE:** Accora Ltd, 38 Main Street, Swords, Co. Dublin, Ireland, K67 E0A2  
**T:** +353 (0)1 695 0614 **E:** info@accora.care **W:** accora.care

® and ™ are trademarks of Accora Limited.  
As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.  
All information is correct at time of printing and is subject to change.