

Configura[®]RiserRecliner

by Accora



Assessment form

Size

Chair size	Part number	Seat depth	Seat height	Seat width	Weight capacity
Extra Small	CFXSD161815DT	400 mm (16")* <input type="checkbox"/>	375 mm (15") <input type="checkbox"/>	400 mm (16")	20 stone
		450 mm (18") <input type="checkbox"/>	400 mm (16")* <input type="checkbox"/>		
			425 mm (17")** <input type="checkbox"/>		
Small Tall	CFSTD182018DT	450 mm (18")* <input type="checkbox"/>	450 mm (18") <input type="checkbox"/>	450 mm (18")	25 stone
		500 mm (20") <input type="checkbox"/>	475 mm (19")* <input type="checkbox"/>		
			500 mm (20")** <input type="checkbox"/>		
Medium Low	CFMLP202016DT	450 mm (18")* <input type="checkbox"/>	400 mm (16") <input type="checkbox"/>	500 mm (20")	20 stone
		500 mm (20") <input type="checkbox"/>	425 mm (17")* <input type="checkbox"/>		
			450 mm (18")** <input type="checkbox"/>		
Large Low	CFLLP222016DT	450 mm (18")* <input type="checkbox"/>	400 mm (16") <input type="checkbox"/>	550 mm (22")	20 stone
		500 mm (20") <input type="checkbox"/>	425 mm (17")* <input type="checkbox"/>		
			450 mm (18")** <input type="checkbox"/>		
Extra Large	CFXLD242018DT	450 mm (18")* <input type="checkbox"/>	450 mm (18") <input type="checkbox"/>	600 mm (24")	25 stone
		500 mm (20") <input type="checkbox"/>	475 mm (19")* <input type="checkbox"/>		
			500 mm (20")** <input type="checkbox"/>		
Extra Large Low	CFXLLP242016DT	450 mm (18")* <input type="checkbox"/>	400 mm (16") <input type="checkbox"/>	600 mm (24")	20 stone
		500 mm (20") <input type="checkbox"/>	425 mm (17")* <input type="checkbox"/>		
			450 mm (18")** <input type="checkbox"/>		

Pressure



Visco Cushion (std) ☐



Castellated ☐



Cushionair ☐



Integrated self supply ☐

Posture



Multi-adjustable Pillow Backrest (std) ☐



2 x Lateral Supports ☐



Lateral Support Backrest ☐



Adjustable Lateral Support Backrest ☐



Cut-down Lateral Support System ☐

Extras

Armrests	Select	Rise setting	Select	Hand control	Select
Overhang on outside (std)	<input type="checkbox"/>	Rise and tilt (std)	<input type="checkbox"/>	RHS (std)	<input type="checkbox"/>
Overhang on inside - reduces width between the armrests by 35 mm (1.5")	<input type="checkbox"/>	Medium rise	<input type="checkbox"/>	LHS	<input type="checkbox"/>
		Flat rise	<input type="checkbox"/>	Carer position	<input type="checkbox"/>

Note - Medium Rise and Flat Rise are only available on Medium and Large chairs not Small and Low Profile chairs

Chair assessment and delivery environmental questionnaire

When carrying out an assessment please fill out the below form prior to ordering the chair.
Please give as much detail as possible.

Access

Is there parking for a large van? yes ☐ no ☐

Are there steps to the property? yes ☐ no ☐

Is there easy access to the property? yes ☐ no ☐

Infection control

Are there pets? yes ☐ no ☐

Is the house smoke free? yes ☐ no ☐

Internal

Are the doors and corridors wide enough for a chair to go through? yes ☐ no ☐
Min. width for most chairs is 650 mm.

Does the chair have to go upstairs? If yes, please specify number of steps and if a stair lift is fitted.
.....

Is space being cleared for the chair? yes ☐ no ☐

Is there a suitable power supply for the chair? yes ☐ no ☐

Other

Does the driver need to be aware of any other points? yes ☐ no ☐
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....