

Configura® Lite 2

The entry-level chair with tilt-in-space

The Lite 2 is the most cost-effective chair in the **Configura** range, featuring a single motor tilt-in-space mechanism, height and depth adjustment and a high-risk pressure cushion as standard.



Low



Key features:



Standard rise



Tilt-in-space with legrest elevation



Adjustable seat height



Adjustable seat depth

Chair size	Part number	Seat depth	Seat height	Seat width	Weight capacity
Small	CHAIR-0-CM2-140	450 mm (18")	400 mm (16")	450 mm (18")	25 stone
		500 mm (20")	450 mm (18")**		
Medium	CHAIR-0-CM2-240	450 mm (18")	450 mm (18")	500 mm (20")	25 stone
		500 mm (20")	500 mm (20")**		
Large	CHAIR-0-CM2-340	450 mm (18")	450 mm (18")	550mm (22")	25 stone
		500 mm (20")	500 mm (20")**		